

# NEW PRODUCT

## Splash

# WHEATON CryoFile™ Freezer Boxes

- Grid ID numbered from 1 – 81 for easy content identification
- The lid and the bottom of the box is numbered and coordinated with the sample ID
- 6 different colors provide preference and identification of the samples
- Water resistant CryoFiles allow for longer durability
- Use with 1, 2 and 5 mL cryogenic vials

Wheaton CryoFile™ and CryoFile XL™ cryogenic vial freezer boxes offer secure storage for high value samples when immersed at 196°C.

The NEW design provides compartments numbered from 1-81 which aids in identifying the contents of the box without the need to check the vials. The numbered grid on the inside lid provides a space for written identification. The bottom of the box uses the same numbering system. Partitions are removable for storage of containers other than cryogenic vials such as 10 mL and 50 mL conical tubes.

The water resistant freezer boxes come in six colors (blue, green, pink, red, white, and yellow) and can withstand cryogenic temperatures up to -196°C. Available in two sizes, the CryoFiles store 1- 2 mL vials, the CryoFiles XL store 5 mL vials vertically and 10 mL and 50 mL horizontally.

Using CryoFiles allow the cryogenic samples to be colored coordinated with the Wheaton Cryule® Cryogenic Vials, which come in the same colors.

See the Wheaton On-Line Catalog for a complete listing of Wheaton cryogenic vials. See the Wheaton On-Line Catalog for a complete listing of Wheaton cryogenic vials.

## Ordering Information

Cat. No.	Description	For Vial Size (mL)	Color	Qty/Cs
<b>CryoFile™ Freezer Box</b>				
W651600	CryoFile Freezer Box	2	Green	15
W651601	CryoFile Freezer Box	2	Yellow	15
W651602	CryoFile Freezer Box	2	Pink	15
W651603	CryoFile Freezer Box	2	White	15
W651604	CryoFile Freezer Box	2	Blue	15
W651605	CryoFile Freezer Box	2	Red	15

## CryoFile XL™ Freezer Box

W651600-XL	CryoFile XL Freezer Box	5	Green	15
W651601-XL	CryoFile XL Freezer Box	5	Yellow	15
W651602-XL	CryoFile XL Freezer Box	5	Pink	15
W651603-XL	CryoFile XL Freezer Box	5	White	15
W651604-XL	CryoFile XL Freezer Box	5	Blue	15
W651605-XL	CryoFile XL Freezer Box	5	Red	15

CryoFile and CryoFile XL are trademarks of Wheaton Industries Inc. Stylized "W" is a registered trademark of Wheaton Industries Inc.

©Wheaton Industries Inc. 2008



## Box Dimensions:

CryoFile:..... 5 x 5 x 1 ¾ in (127 x 127 x 44 mm)

CryoFile XL:..... 5 x 5 x 3 ¾ in (127 x 127 x 95 mm)

## Removable Partition Dimensions:

CryoFile:.....4 7/8 x 4 7/8 x 1 1/8 in  
(123.8 x 123.8 x 28.6 mm)

CryoFile XL:.....4 7/8 x 4 7/8 x 2 1/8 in  
(123.8 x 123.8 x 60.3 mm)



**WHEATON**  
Science International  
a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038 USA  
+1.856.825.1100 ■ +1.856.825.1368 (F)

[www.wheatonsci.com](http://www.wheatonsci.com)