

FOR FURTHER INFORMATION,
CONTACT JACKIE WILLIAMS
+1.856.825.1100, ext. 2889
Jackie.Williams@wheaton.com

FOR IMMEDIATE RELEASE
DATE: MARCH 7, 2009

**REDUCE HAND FATIGUE AND IMPROVE PRODUCTIVITY WITH THE NEW WHEATON
CRIMPENSTEIN ELECTRIC CRIMPER/DECAPPER**

Millville, NJ – Wheaton Science Products introduces the new Crimpenstein Electric Crimper/Decapper. This lightweight, semi-automatic product is ideal for laboratories that crimp or decap large volumes of aluminum seals. Its ergonomic design also helps to reduce hand fatigue caused by the use of manual crimpers and decappers.

The Crimpenstein is available with a choice of four different size crimper or decapper jaws (8, 11, 13 and 20 mm). The ability to quickly change between different sizes or between crimping and decapping offers an unmatched level of flexibility and economy for applications where various sizes of aluminum seals are required. The Crimpenstein features an adjustable jaw intensity which allows the user to regulate jaw action as required to achieve the perfect crimp. A memory feature also retains the last crimp setting when changing from decapping back to crimping. Since it is AC powered (electric), there is no need for regulated compressed air or batteries which discharge over time. A sturdy stand is available which will hold the Crimpenstein in a fixed vertical position for operations where vials need to be handled one at a time.

Wheaton Science Products is a division of Wheaton Industries Inc. Wheaton is a fully integrated business offering a large product portfolio of commodity and leading edge products and services for the pharmaceutical and life science R&D, diagnostic packaging, and specialty chemicals market. Our products are used in the up-front stages of research and development and later in full scale production.

For more information, please visit our website www.wheatonsci.com or contact Wheaton Science Products Customer or Technical Service at 800.225.1437 or +1.856.825.1100, Fax: +1.856.825.9035, e-mail: Jackie.Williams@wheaton.com or write 1501 North 10th Street, Millville, NJ 08332-2038.

###