

FOR FURTHER INFORMATION,  
CONTACT JACKIE WILLIAMS  
+1.856.825.1100, ext. 2889  
Jackie.Williams@wheaton.com

**FOR IMMEDIATE RELEASE**  
DATE: MARCH 7, 2009

## **NEW Wheaton CryoFile™ and CryoFile XL™ for Cryogenic Storage**

Millville, NJ – Wheaton Science Products announces the introduction of the CryoFile and CryoFile XL cryogenic storage boxes. The NEW Wheaton CryoFile and CryoFile XL will save researchers' time identifying samples and offer secure cryogenic storage for samples to -196°C.

The new design provides an easy ID grid numbered from 1-81 which aids in identifying the contents of the box without the need to check the vials. The numbered grid on the inside lid provides a space for written identification information. The bottom of the box has the same numbering system. The partitions in the box are removable for storing contents other than cryogenic vials such as 10 mL and 50 mL conical tubes horizontally.

The water resistant freezer boxes come in six colors (blue, green, pink, red, white, and yellow) Wheaton Cryule® Cryogenic Vials are available in the same colors for easy coordination with the correct freezer box. The CryoFiles are available in two sizes; the CryoFiles store 1- 2 mL vials and the CryoFiles XL store 5 mL vials vertically and 10 mL and 50 mL horizontally.

Wheaton Science Products is a division of Wheaton Industries Inc. Wheaton is a fully integrated business offering a large product portfolio of laboratory glass and plastic ware and leading edge products and services for the pharmaceutical and life science R&D, diagnostic packaging, and specialty chemicals market. Our products are used in the up-front stages of research and development and later in full scale production.

For more information, please visit our website [www.wheatonsci.com](http://www.wheatonsci.com) or contact Wheaton Science Products Customer or Technical Service at 800.225.1437 or +1.856.825.1100, Fax: +1.856.825.9035, e-mail: [Jackie.Williams@wheaton.com](mailto:Jackie.Williams@wheaton.com) or write 1501 North 10<sup>th</sup> Street, Millville, NJ 08332-2038.

###