



Company Contact: Tracy Neri-Luciano
WHEATON
+1.856.825.1100
Email: tracy.neri-luciano@WHEATON.com
www.WHEATON.com

Agency Contact: Chris Hoag
Kenyon Hoag Associates
+1.201.236.9898
Email: choag@kenyonhoag.com
www.kenyonhoag.com

WHEATON introduces the new SingleScan™ Bar Code Reader

Sample traceability and enhanced security is facilitated with reliable, plug-and-play technology

Adding to our portfolio of new and innovative products developed specifically for biological sample management, WHEATON is pleased to introduce the new SingleScan Bar Code Reader. This is the latest introduction for WHEATON, advancing our charter to design and deliver today's most innovative solutions for the laboratory research, diagnostic packaging and specialty pharmaceuticals industries.

The WHEATON SingleScan Bar Code Reader was created as a tool to aid researchers in making sample traceability more efficient and accurate. A sophisticated 2D Data Matrix and 1D Linear Bar Code Reader capable of reading coded vials and ampules from a wide range of manufacturers and label printers, the new SingleScan Bar Code Reader offers enhanced reading performance in a compact unit. Researchers simply hold or place the vials, such as our CryoELITE® Cryogenic vials, over the scanning element for rapid reading and precise identification.

The SingleScan features a user-friendly plug-and-play design that enables an easy set-up with no software installation required; simply plug into the USB port and open the application file. The USB interface allows for an easy connection to any computer and enables multiple user/multiple software applications. Additionally, the inputs decoded bar code ID into any software application where the cursor is placed, making the SingleScan the easiest and most cost effective way to manage and track biological samples.

"We are seeing more and more interest in productivity in laboratories, especially those labs performing biobanking," says Stephen R. Drozdow, President of WHEATON. "Our approach has always been to make the researcher's day-to-day tasks easier and more effective letting them concentrate on their life's work," he added.

The SingleScan Bar Code Reader is another addition to WHEATON BioBanking Solutions™, a comprehensive offering of high quality products and customized services for research that also includes glass or plastic cryogenic containers, biobank storage boxes, and barcoding and scanners. For more information, please visit <http://www.WHEATON.com> or contact WHEATON at 800.225.1437 in the USA, or 856.825.1100 for international callers.

About WHEATON

WHEATON admires those people who devote their life to science. It's this admiration that drives us to design and deliver today's most innovative solutions for the laboratory research, diagnostic packaging and specialty pharmaceuticals industries. It's this admiration that motivates us to offer the best customer service experience and most dependable products in the industry. And, it's with admiration that we stand behind those people who help move the world forward one discovery, one sample at a time. The WHEATON name means more than a legacy of the highest quality products and services; it stands for a commitment to honor the effort of all scientists... because it's their life's work.