


Narrow Mouth



### Media Lab Bottles, Graduated

- Clear borosilicate glass bottle with graduations and marking patch
- Amber borosilicate glass bottle in 125mL size only
- Pre-attached black phenolic screw caps with choice of three liners
- Hand grips on 1L bottle provide easy handling

Cat. No.	Capacity (mL)	Cap On Dia x H (mm)	Cap Off Dia x H (mm)	Cap Size	Qty/Case
<b>Clear Glass Bottle Without Cap</b>					
219435	125	55 x 123	55 x 119	33-430	48
219437	250	67 x 152	67 x 148	33-430	48
219439	500	88 x 192	88 x 188	33-430	24
219440	1000	103 x 235	103 x 231	38-430	24
<b>Clear Glass Bottle With Polyethylene (LDPE) Lined Phenolic Cap</b>					
219715	125	55 x 123	55 x 119	33-430	48
219717	250	67 x 152	67 x 148	33-430	48
219719	500	88 x 192	88 x 188	33-430	24
219720	1000	103 x 235	103 x 231	38-430	24
<b>Clear Glass Bottle With Rubber Lined Phenolic Cap (Autoclavable)</b>					
219755	125	55 x 123	55 x 119	33-430	48
219757	250	67 x 152	67 x 148	33-430	48
219759	500	88 x 192	88 x 188	33-430	24
219760	1000	103 x 235	103 x 231	38-430	24
<b>Clear Glass Bottle With PTFE Faced Rubber Lined Phenolic Cap (Autoclavable)</b>					
219815	125	55 x 123	55 x 119	33-430	48
219817	250	67 x 152	67 x 148	33-430	48
219819	500	88 x 192	88 x 188	33-430	24
219820	1000	103 x 235	103 x 231	38-430	24
<b>Amber Glass Bottle Without Cap</b>					
219845	125	55 x 123	55 x 119	33-430	48
<b>Amber Glass Bottle With Polyethylene (LDPE) Lined Phenolic Cap</b>					
219855	125	55 x 123	55 x 119	33-430	48
<b>Amber Glass Bottle With Rubber Lined Phenolic Cap (Autoclavable)</b>					
219865	125	55 x 123	55 x 119	33-430	48
<b>Amber Glass Bottle With PTFE Faced Rubber Lined Phenolic Cap (Autoclavable)</b>					
219875	125	55 x 123	55 x 119	33-430	48



### Media Lab Bottles, Non-Graduated

- Clear borosilicate glass bottle without graduations
- Amber borosilicate glass bottle in 125mL size only
- Pre-attached black phenolic screw caps with choice of liner
- Hand grips on 1 Liter bottle provide easy handling

Cat. No.	Capacity (mL)	Cap On Dia x H (mm)	Cap Off Dia x H (mm)	Cap Size	Qty/Case
<b>Clear Glass Bottle Without Cap</b>					
219415	125	55 x 123	55 x 119	33-430	48
219417	250	67 x 152	67 x 148	33-430	48
219419	500	88 x 192	88 x 188	33-430	24
219420	1000	103 x 235	103 x 231	38-430	24
<b>Clear Glass Bottle With Polyethylene (LDPE) Lined Phenolic Cap</b>					
219495	125	55 x 123	55 x 119	33-430	48
219497	250	67 x 152	67 x 148	33-430	48
219499	500	88 x 192	88 x 188	33-430	24
219500	1000	103 x 235	103 x 231	38-430	24
<b>Clear Glass Bottle With (Autoclavable) Rubber Lined Phenolic Cap</b>					
219575	125	55 x 123	55 x 119	33-430	48
219577	250	67 x 152	67 x 148	33-430	48
219579	500	88 x 192	88 x 188	33-430	24
219580	1000	103 x 235	103 x 231	38-430	24
<b>Amber Glass Bottle without Cap</b>					
219885	125	55 x 123	55 x 119	33-430	48

### Replacement Caps for Media Bottles

- For WHEATON media bottles
- Black phenolic caps with 33 and 38 mm screw threads

Cat. No.	Description	Size (mm)	Qty/Case
240080	Solid Cap with Polyethylene (LDPE) Liner	33-430	200
240081	Solid Cap with Polyethylene (LDPE) Liner	38-430	200
240280	Solid Cap with Rubber Liner (Autoclavable)	33-430	200
240281	Solid Cap with Rubber Liner (Autoclavable)	38-430	200
240480	Solid Cap with PTFE Faced Rubber Liner (Autoclavable)	33-430	100
240481	Solid Cap with PTFE Faced Rubber Liner (Autoclavable)	38-430	100
240680	Open Top Cap with Gray Chlorobutyl/50 Septa and Flange (Autoclavable)	33-430	100
240683	Open Top Cap with Gray Chlorobutyl/50 Septa and Flange (Autoclavable)	33-430	1000



Solid Top Cap with Liner



Open Top Cap with Septa

# WHEATON

## Wide Mouth Media Bottles



### Media Bottle, Lab 45™ Graduated Bottles

- Clear Borosilicate glass with graduations
- Available in 100mL - 2L size with wide mouth
- Graduated with writing patch
- Special pour lip facilitates pouring
- Provided with or without screw caps
- Color coded caps can be purchased separately

Cat. No.	Capacity (mL)	Dia x H (mm)	Cap Size (mm)	Qty/Case
<b>Bottles with Caps and Pour Rings Attached</b>				
219925	100	56 x 109	45	12
219927	250	70 x 144	45	12
219929	500	86 x 186	45	12
219930	1000	101 x 233	45	12
219931	2000	137 x 267	45	6
<b>Bottles without Caps or Pour Rings</b>				
219915	100	56 x 103	45	12
219917	250	70 x 138	45	12
219919	500	86 x 180	45	12
219920	1000	101 X 227	45	12
219921	2000	137 X 261	45	6



### Replacement Caps for Lab 45™ Media Bottles

- 45mm screw cap with or without sealing ring
- High temperature screw caps
- Caps and stoppers for lyophilization
- Pour rings

The 45mm polypropylene screw cap with sealing ring are replacements for caps that are provided with Lab 45 media bottles. For high temperature applications use the PBT (polybutylene terephthalate) screw caps and the ETFE (ethylene tetrafluoroethylene) pour lips which can withstand temperatures up to 180°C. For lyophilization applications use the solid or open top polypropylene screw cap without the sealing ring and use the black lyophilization stopper.

Cat. No.	Description	Size (mm)	Qty/Case
240726	White PP cap with sealing ring	45	12
240726-03	Red PP cap with sealing ring	45	12
240726-04	Blue PP cap with sealing ring	45	12
240736	White PP cap without sealing ring	45	12
240740	White PP cap with PTFE faced silicone liner	45	12
240746	White PP cap with open top	45	12
240750	Red PBT cap with PTFE Faced Silicone Liner	45	10
240756	Natural Polypropylene Pour Rings	45	12
240760	Red ETFE Pour Rings	45	10
224100-503	Halobutyl/60 Black Lyophilization Stoppers	43	12
224100-507	EPDM Black Lyophilization Stoppers	43	12

Lab 45 is a trademark of WHEATON Industries Inc. The stylized W is a registered trademark of WHEATON Industries.



### Media Bottle, Sterile PET

- New improved 1 Liter Bottle
- Thicker wall weight
- Wide mouth design
- Permanent in-mold graduations
- Linerless closure
- Caps attached

Cat. No.	Size (mL)	Dimensions (L x W x H)	Cap Size	Qty/Tray	Qty/Case
W219980	1000	91 x 91 x 220 mm	38-430	12	24



125 & 500 mL Bottles

### Media Bottle, Sterile PET

- Permanent in-mold graduations
- No-drip pour lip
- Lightweight
- 20% headspace for additives

WHEATON PET Square Media Bottles are manufactured from polyethylene terephthalate which provides the clarity of glass, the strength of plastic, and a space-saving design. These sterile bottles are tested superior for pH stability, temperature durability, cloning efficiency and cytotoxicity.

Cat. No.	Size (mL)	Dimensions (L x W x H)	Cap Size	Qty/Tray	Qty/Case
219975	125	53 x 53 x 99 mm	33-430	24	48
219979	500	75 x 75 x 175 mm	33-430	12	24



**WHEATON Science Products**

a division of WHEATON Industries Inc.

1501 North 10th Street, Millville, NJ 08332-2038

800.225.1437 (U.S. & Canada) • 856.825.1100 • 856.825.1368 (F)

[www.wheatonsci.com](http://www.wheatonsci.com)