

Plug Styles

The following Plug Styles can be used with Wheaton magnetic stirrers, roller culture apparatus, incubators and pumps.



Plug Code "A": North America Plug



Bahamas	Guatemala	Philippines
Bolivia	Haiti	Puerto Rico
Brazil	Honduras	Saudi Arabia
Canada	Jamaica	Singapore
Costa Rica	Korea, S	Sweden
Dominican Rep.	Libya	Taiwan
El Salvador	Nicaragua	United States
Ecuador	Norway	Venezuela
Guam	Panama	Virgin Islands
Peru		




Plug Code "B": Japan Plug




Bolivia
Japan




Plug Code "C": Continental Europe Plug



Austria	Indonesia	Singapore
Belgium	Latvia	Sweden
Bolivia	Lithuania	Spain
Bulgaria	Luxembourg	Tahiti
Czech. Rep	Mozambique	Turkey
Estonia	Netherlands	USSR (former)
Finland	Poland	Vietnam
France	Portugal	Yugoslavia (former)
Germany	Romania	Zaire (Rep. of)
Greece	Russia	
Hungary	Slovakia	
Iceland	Slovenia	




Plug Code "D": United Kingdom Plug



Bahrain	Kuwait	Tanzania
Belgium	Malaysia	Trinidad
Bermuda	Mauritius	& Tobago
Cyprus	Malta	United Arab
England	Nigeria	Emirates
Myanmar	N. Ireland	United Kingdom
Ghana	Oman	Wales
Hong Kong	Qatar	Zambia
Ireland	Scotland	Zimbabwe
Iraq	Singapore	Zimbabwe
Jordan	Sudan	
(Rep. of)		



Plug Code "F": Australia/China Plug




Argentina
Australia
China (People's Rep. of)
Fiji
New Zealand
Papua New Guinea
Uruguay



Plug Code "G": Italy/Chile Plug



Chile
Italy




Plug Code "J": India Plug

Bangladesh
India
Pakistan
South Africa
Sri Lanka
Uganda

Special Plug Styles




Plug Code "I": Denmark Plug



Denmark
Greenland



Plug Code "K": Israel Plug



Israel
Jamaica



Plug Code "L": Switzerland Plug



Switzerland



Additional plug styles are available through special order. Please contact Technical Service for additional information.