

WHEATON

**Purchase 2 Cases of
WHEATON Liquid Scintillation
Vials and Receive 1 Case FREE***



Liquid Scintillation Vials

- 20mL size
- Made from WHEATON 180 low potassium borosilicate glass
- Background counts are consistent and low; ultraviolet transmission is high
- Screw caps attached to vials

Cat. No.	Cap Size	Cap Material	Liner Material	Qty/Case
986540	22-400	Polypropylene	Foamed Polyethylene	500
986541	22-400	Polypropylene	Metal Foil	500
986542	22-400	Urea	Metal Foil	500
986546	22-400	Urea	Polyethylene Cone	500
986548	22-400	Urea	Polyethylene Disc	500
986560	24-400	Polypropylene	Foamed Polyethylene	500
986561	24-400	Polypropylene	Metal Foil	500
986562	24-400	Urea	Metal Foil	500
986568	24-400	Urea	Polyethylene Disc	500

**Complete the redemption coupon below
and send or mail/fax/email it along with
your proof of purchase to:**

Limited Time Offer
Offer Valid thru June 30, 2012

WHEATON

Attention: Joel Torres
1501 North 10th Street
Millville, NJ 08332-2038 USA

Email: joel.torres@wheaton.com
800.225.1437 (T) | 856.825.9035 (F)

Please send my FREE case of WHEATON Liquid Scintillation Vials Cat. No. _____.

Please enclose or fax a copy of your proof of purchase: Invoice, packing slip, or purchase order with your redemption.

Name: _____
Company: _____
Address: _____
Address: _____
City: _____
State: _____
Zip/Postal: _____
Country: _____
Tel: _____
E-Mail: _____

This offer is void where prohibited by law or company policy. All information must be completed to be a valid entry. Offer valid in North America only thru 06/30/2012. Limited to two redemptions per customer. *Free case must be of equal or lesser value. Not for resale.

www.wheaton.com